







DO NOT LEAN

nt
ultraschall

steinberger medizintechnik

**Siemens Acuson
X300PE
#315988 SMT#3276**



SIEMENS

Siemens Medical Electronics AG, Postfach 1015619, 69126 Heidelberg, Germany
Siemens AG, Postfach 1015619, 69126 Heidelberg, Germany

Siemens AG, Postfach 1015619, 69126 Heidelberg, Germany
Siemens AG, Postfach 1015619, 69126 Heidelberg, Germany

Note: See operating instructions.

Hinweis: Gebrauchsanweisung beachten.

Remarque: Se reporter au Mode d'emploi.

Nota: Ver instrucciones de uso.

Nota: Consultare il manuale dell'utente.

注意: 取扱い説明書を読むこと

주의: 사용자 설명서를 반드시 읽어 주십시오.

注意: 请参照使用说明书。


Caution: Federal (U.S.) law restricts this device to sale by
on the order of a physician.

Diagnostic Ultrasound System

ACUSON X300 PREMIUM EDITION

Model 10348531

Manufacture Date Jan. 2010



(1P)Model No. 10348531




(S)Serial No. 315988


Product Status Number 02

Siemens Medical Solutions USA, Inc. 1230 Shorebird Way.

Mountain View, CA94043 USA.


Made in Korea

 System
IVK



Upgrade

A B C D E F G H J K L M N P R S T V W



Steckbrief:

Siemens Acuson X300 PE Farbsystem

BJ Jan. 2010

18" LCD-Monitor

3 Sondenports

digitaler Bildspeicher auf Festplatte

vordefinierte Messprogramme

PW-Doppler, Power Doppler, opt. CW-Doppler

anatom. M-Mode optionales Zubehör

Bildverbesserungssoftware:

DTCE™ (Dynamic tissue contrast enhancement Technology)

SieClear™ (Multi-View-spatial compounding)

THI (Tissue Harmonic Imaging)

Multi-Beam Formation

Interne Nr.: 3276

Preis 5200€ incl. MwSt. (Endpreis) incl. Abdomensonde

Weitere Sonden auf Anfrage

Optional: Printer

Garantie: 6 Monate

Generalüberholt

Service und Wartung neu

Systemstandort: München